



BlueClassic Prestained Protein Marker

10-180 kDa
Ready-to-use prestained protein ladder

Cat. No.	Amount
PS-107	2 x 250 µl



Migration Pattern

For in vitro use only!

Shipping: shipped on blue ice

Storage Conditions: store at -20 °C

Short term storage (up to 1 month) at 4 °C possible.

Shelf Life: 12 months

Purity: > 95 % (SDS-PAGE)

Form: liquid (Supplied in 20 mM Tris-phosphate pH 7.5 at 25°C, 2 % SDS, 1 mM β-mercaptoethanol, 3.6 M urea and 15 % glycerol)

Concentration: 1 mg/ml

Description:

The BlueClassic Prestained Protein Marker is a blue-colored protein standard with 10 pre-stained and purified proteins supplied in gel loading buffer for direct application in SDS-polyacrylamide gel.

Recommended Load:

3-5 µl Protein Molecular Weight Marker per lane